



Using the IQM transmission detector to monitor linac performance

- Dan Pavord
- Chief of Medical Physics
- Allegheny Health Network





- Conflict of interest
- I have received a grant from Elekta to investigate the use of IQM for patient specific QA.




The IQM device (Integral Quality Monitor) iRT Systems, Koblenz, Germany



EMBRACING CHANGE.
IMPACTING PATIENT CARE.

AAPM 2024
66TH ANNUAL MEETING & EXHIBITION





Patient Specific QA

- Ghafarian, M., Price, M., & Morales-Paliza, M. (2021). Comparison of pretreatment VMAT quality assurance with the integral quality monitor (IQM) and electronic portal imaging device (EPID). *Journal of applied clinical medical physics*, 22(3), 166-175.
- Saito, M., Sano, N., Shibata, Y., Kuriyama, K., Komiyama, T., Marino, K., ... & Onishi, H. (2018). Comparison of MLC error sensitivity of various commercial devices for VMAT pre-treatment quality assurance. *Journal of Applied Clinical Medical Physics*, 19(3), 87-93. (Delta4, PTW 2D array, IQM, farmer chamber)



Real Time Monitoring of Treatment Delivery

- Casar, B., Pasler, M., Wegener, S., Hoffman, D., Talamonti, C., Qian, J., ... & Peterlin, P. (2017). Influence of the integral quality monitor transmission detector on high energy photon beams: a multi-centre study. *Zeitschrift Für Medizinische Physik*, 27(3), 232-242.
- The effects of the IQM detector on photon beam properties were found to be small yet statistically significant. The magnitudes of changes which were found justify treating IQM either as tray factors within the treatment planning system (TPS) for a particular energy or alternatively as modified outputs for specific beam energy of linear accelerators, which eases the introduction of the IQM into clinical practice.
- Arilli, C., Wandael, Y., Galeotti, C., Marrazzo, L., Calusi, S., Grusio, M., ... & Talamonti, C. (2020). Combined use of a transmission detector and an EPID-based in vivo dose monitoring system in external beam whole breast irradiation: a study with an anthropomorphic female phantom. *Applied Sciences*, 10(21), 7611.
- IQM can detect output and MLC/jaw positioning errors. EPID can detect output and patient positioning errors.



Linac QA

Received: 27 June 2023 | Revised: 9 December 2023 | Accepted: 14 December 2023

DOI: 10.1002/acm2.14245

RADIATION MEASUREMENTS

JOURNAL OF APPLIED CLINICAL
MEDICAL PHYSICS

Feasibility of using a transmission ion chamber for QA tests of medical linear accelerators

Makan Farrokhkish¹ | Andrew J. Veres² | Luis E. Fong de los Santos² |
John J. DeMarco³ | Mohammad K. Islam^{1,4}

Output
Symmetry
Flatness
MLC position
Dose Rate
MU Linearity




Linac QA, photons

- Add in more checks
 - Energy
 - Jaw Position
 - Wedge Factor
 - Gantry/Collimator Angle
 - Beam Center

Output
Symmetry
Flatness
MLC position
Dose Rate
MU Linearity




Automated Field Sequence, Elekta Versa HD, 6MV






Automated Field Sequence






Automated Field Sequence






Automated Field Sequence

Repeat at four cardinal gantry angles





Output Constancy



Note: On Day 13 deliberate +/-1% changes made to output.




Energy

- Adjusted energy as determined by PDD at 10cm
- Measured resulting change in ratio of small/large field signal with IQM to establish Action Level




Energy







Energy






Jaw Position, use wedge to create gradient





Jaw Position, alternate method of Kilian Michel, iRT



- If Y1 was moved by 1mm it would only affect the IQM signal by less than 1%, because the relating radiation contribution would be absorbed by the closed MLC-pair.
- Similarly, the position/opening of the closed MLC-pair has a negligible effect on the IQM signal because it is half-way covered by Y1.
- However, if Y2 moved, it would have a dramatic effect on the open field size and thereby on the IQM Signal. (preliminary results ~11%/mm)



Wedge Factor

Field	IQM WF	measured
10x10	0.289	0.268
20x20	0.308	0.279



Gantry/Collimator Angle

Data easily captured from log files.

All within 0.5 degree of expected angles.





Beam Center

Ratios of $+X/-X$ and $+Y/-Y$ have been very stable ($\pm 0.5\%$).
Next step is to purposely steer beam center and see if ratios correlate.



Electrons





Electrons

- Seven Measurements
 - 20x20 with energy filter
 - 20x20 open
 - 4x4 at CAX
 - Offset field at 4 collimator angles
- Beam analysis
 - 20x20 open for output
 - Ratio 20x20 filter/open for energy
 - Ratios of offset fields to 4x4 CAX for flatness/symmetry




Electrons Energy

Energy Ratio Measured with IQM					
Day	6MeV	8MeV	10MeV	12MeV	15MeV
0	0.405	0.533	0.668	0.777	0.905
7	0.413	0.537	0.674	0.784	0.913
8	0.404	0.540	0.675	0.784	0.911
9	0.405	0.540	0.674	0.782	0.909




Electrons Energy





Electrons Energy





Time

- Currently about 70 minutes for all tests, photon and electron
- Compares well with 90 minutes using our existing methods



- Acknowledgements
- Alex Solodkin, Elekta
- Jurgen Oellig, iRT Systems
- Will Roestel, iRT Systems
- Kilian Michel, iRT Systems